

A reliable
family business,
since 1891

OTOLIFT AIR



reddot design award
winner 2011





CONTENTS

Introduction	4
OTOLIFT AIR	6
Six benefits	12
Rail system	14
Equipment	18
Choice of colours	24
Installation	26
About Otolift	27

Otolift – more than a name,
a promise too!



*'Every day I'm grateful for
choosing Otolift.'*

Mr Dunning, Warwick

Over 150,000 customers already know that an Otolift stairlift is much more than just an accessory in people's daily lives. The models manufactured by the Dutch market leader - Otolift; the first of which was produced back in 1968, enable people to continue living life in the manner they are accustomed to, self-supporting and independent, with pleasure and joy in life.



Otolift customers would like to spend their lives in familiar surroundings, the place they have called home for years, or even decades. Long time friends, neighbours and relatives are close by and they can still carry out the routines and hobbies locally that make life worthwhile.

These are things that people can continue to enjoy even if a staircase restricts their mobility.

Choosing an Otolift Stairlift assures you of not only excellent service and outstanding product quality, but also allows you to live the life you want to lead.

Our customers confirm this year in, year out; Otolift offers a degree of freedom. The freedom that they deserve.

OTOLIFT AIR

True innovation



Following intense research and extensive materials tests, we have succeeded in developing a trendsetting rail: the 6cm diameter rail is currently* available exclusively for the Otolift Air. The Otolift Air has become our customers' favourite, thanks to its modern design and advanced technology.



This ground-breaking technology ensures that you experience Otolift's proven convenience and its accredited safety, while also allowing you to revel in its surprisingly modest elegance. A single, unobtrusive rail fitted close to the stairs ensures the complete rail blends harmoniously with your living environment.



The complete installation evokes enthusiastic reactions - even if they don't notice it at first glance.



reddot design award
winner 2011

*01-03-16

OTOLIFT AIR

New elegance

The Otolift rail requires very little space on the staircase.

Mrs Dunning, Warwick



The Otolift Air impresses with its modern design, which does not immediately capture one's attention, but merges into the surroundings in an elegant and modest manner.



The Otolift Air evokes surprise when viewed more closely: while its unique rail appears rather delicate, it easily surpasses the Otolift ONE in terms of possibilities and durability. Elegance and technology are perfectly combined in this system.

The 6cm diameter rail is available exclusively from Otolift Dealers nationwide.



The Otolift Air impresses people with its exclusive elegance, excellent comfort and vital safety functions.

OTOLIFT AIR

Flexible convenience

Retired Architect, Paul Taylor, of Guildford, still has a keen eye for detail. He therefore kept a close watch during the installation of his Otolift Air.

“I was very much impressed with the installation of my stairlift. Two engineers took just three hours to complete the install without any modifications to my stairs. They left my house spotlessly clean and gave a detailed handover, ensuring I could operate the lift safely before they left.”



‘The Otolift Air does not hinder the use of the staircase in any way.’

Mr Taylor, Guildford

The Otolift Air fits neatly around the inside or outside of your staircase. Fitting it on the inside of the stairs often means your handrail can remain on the opposite wall, leaving the widest part of the stairs clear for others to use safely.



Steep, narrow or spiral staircases. The Otolift effortlessly accommodates the most demanding of staircases and discreetly follows all bends. What's more, the rail requires so little space that normal use of the staircase is not hindered. This ensures not only you, but also your family members or guests can continue to use the staircase safely.

6 benefits of OTOLIFT AIR

- 
1. Your handrail can often remain in place
 2. The widest part of the stairs remains accessible
 3. Unique 6cm diameter rail
 4. Extremely compact chair
 5. Aesthetically pleasing features both folded and unfolded
 6. Powered folding footrest as standard



OTOLIFT AIR

The single-rail system of the future



Choosing the Otolift Air will bring you a product at the forefront of rail design. The 6cm diameter tube used to manufacture the rail is the smallest diameter of all single tube stair-lift rails on the current market*.



*O1-O3-16

End View of Otolift ONE Rail Section

Ø 8 cm



A true 1:1 scale representation of both the Otolift ONE & Otolift AIR rail system. This clearly shows the difference between the smaller, 6cm diameter tube of the AIR over the larger 8cm diameter tube on the Otolift ONE.

Actual size

End View of Otolift AIR Rail Section

Ø 6 cm



Actual size

OTOLIFT AIR

The spatial miracle

*'Our Otolift is so easy to operate
that I use it to go up and down
whenever I choose.'*

Mrs Dunning, Warwick



The Otolift Air is a prime example of advanced technology and modern design, which subtly blends with the latest living trends.

If the Otolift does happen to get in the way, then the chair can be concealed from view with a single click on the remote control. The chair can also be folded in just two easy movements, in order to save more space.

Automatically rotating chair

The Otolift Air's chair can be rotated conveniently and automatically. This enables you to select the most comfortable position when mounting and dismounting.



Drop Nose rail start available offering a minimum intrusion at the bottom of the stairs.



Otolift stairlifts - Safe and reliable

Safe



The safety belt is easy to use and can be secured using just one hand.

Practical



The safety belt is stored neatly when not in use and offers a sense of security when the lift is in use.

User friendly



The hand control is very robust and break-resistant.

Space saving



When not in use your Otolift can be effortlessly folded away.

Intelligent



Your Otolift stops whenever it encounters an obstruction: Automatically, gently and safely.

Convenient



The remote control can be used to either call the lift to you or send it away when not in use.

Regardless of which model of Otolift stairlift you choose, you will be reassured by the range of advantages that only we have to offer. By opting for an Otolift you can rely on an extensively tried and tested stairlift of premium quality, which offers the optimum in convenience and safety, regardless of the circumstances, every single time.

Unique: Fully automatic chair

Flexible



The chair can be rotated automatically to suit even the most challenging staircase...



... giving the flexibility to travel facing in either direction to avoid obstacles.

Comfortable



The footrest can be retracted automatically with just one click on the button.

OTOLIFT AIR

Ideal for steep and narrow staircases



The Otolift Air has no problems with a steep staircase. The seat can also travel in a rotated position to negotiate any staircase, ensuring a more comfortable ride with highest possible clearance for your knees.



Your OTOLIFT AIR Materials of your choice


Not only does a stairlift stand for your independence, but it can also be finished to suit your personal taste and home interior. Modest, or rather more extravagant? Thanks to the various finishes available – including light, dark and natural wood patterns, or even brushed aluminium – the Otolift Air can be tastefully adapted to your interior decor. You also have a choice of upholstery in various high quality fabrics and leather.



Dark wood pattern




Beige
Fabric



Cream
Fabric



Chestnut
Leather




Light wood pattern



Beige
Fabric



Grey
Fabric




Chestnut
Leather




Natural wood pattern



Beige
Fabric



Chestnut
Leather




Brushed aluminium



Taupe
Leather




Chestnut
Leather



Rail colour

Cream
RAL 9001



Light Grey
RAL 7035



Pale Brown
RAL 8025



The actual colours may vary, depending on the print.

OTOLIFT AIR

Installation



→ Neatly completed within 3 hours



→ We only leave after providing calm and clear instructions

→ Screwed to the staircase


About Otolift.

Otolift is a Dutch family company, established in 1891.

Since 1968, Otolift has been 100% producer of top quality stairlifts, and takes a lot of pride in being the benchmark in the industry.



→ 100% Produced by Otolift



27

→ Since 1891

→ 100% Designed by Otolift

Your Otolift approved dealer:

OTOLIFT

WWW.OTOLIFTSTAIRLIFTS.CO.UK
WWW.OTOLIFTSTAIRLIFTS.IE