




OTOLIFT TWO
YOUR DISCREET STAIRLIFT





CONTENTS

4 OTOLIFT TWO The first thing you notice about the TWO, is how little space it occupies. The stairlift is compact but very soundly constructed.	6 YOU CAN STILL USE THE STAIRS When designing the TWO, Otolift did everything possible to make the rail and the chair as compact as possible.	Otolift TWO 4
8 SUPER-SLIM DESIGN The TWO is distinctive for its very slim tubes.	10 NEGOTIATES THE BENDS WITH EASE Are your stairs narrow? Then you will be pleasantly surprised how well the Otolift TWO runs along the curve.	You can still use the stairs 6 Super-slim design 8 Negotiates the bends with ease 10
14 THE CLEVER SOLUTIONS OF THE TWO The TWO offers a good range of clever solutions.	23 ABOUT OTOLIFT Otolift is a Dutch family company, established in 1891.	Intuitive and easy to use 12 The clever solutions of the TWO 14 Choose your upholstery and paint colour 16 Specials 20 Installation 22 About Otolift 23



**"I THOUGHT ALL STAIRLIFTS
WERE ALIKE"**

"What I had not expected was that there was such a difference between stairlifts! Eventually, I chose an Otolift. Why? Because it soon became clear that Otolift was the only company that really paid attention to all details of the lift, up until the supports on the stairs. And I always think that if a company takes so much care about its products, it must be a good company. Also, they had the most beautiful stairlifts."

Mr Harrison, Dublin

OTOLIFT TWO

The first thing you notice about the TWO, is how little space it occupies. The stairlift is compact but very soundly constructed. The Otolift TWO is definitely a stairlift with charm. The subtle, timeless design combines extremely well with any interior.



Advantages of the TWO

- Maximises space on the stairs for non stairlift users
- Strong, safe and reliable build quality
- Blends in with any interior
- Tight bends are possible, thereby minimising loss of space



OTOLIFT **TWO**

You can still use the stairs



When designing the TWO, Otolift did everything possible to make the rail and the chair as compact as possible. The results are clear: The tubes are slim and they can be mounted very close to the wall. The chair is as compact as it is comfortable. As a result, everyone is able to use your stairs with ease.



OTOLIFT **TWO**

Super-slim design



The TWO is distinctive for its very slim tubes. Often these are mounted very close to the wall or banister so that the rail blends in beautifully with the interior and is hardly noticeable at all.



OTOLIFT **TWO**

Negotiates the
bends with ease



Are your stairs narrow? Then you will be pleasantly surprised how well the Otolift TWO runs along the curve. Even at the sharpest bends, the TWO ascends and descends in a single flowing motion. The construction is also extremely sound. So that going up or going down, you always feel perfectly safe.



OTOLIFT TWO

Intuitive and easy to use

STEP 1

Bring the chair towards you using the wireless remote control.



STEP 2

Simply unfold the chair.




STEP 3

Sit down comfortably.



→ SAFE

The chairlift will continue to operate even if there is power cut.



→ CONVENIENT

You can fold the footrest up or down using the conveniently positioned lever.



→ SAFE

The safety belt is easily fastened with one hand.



Your Otolift has been designed for ease of use.
Using it is very intuitive, and will soon become second nature.

STEP 4

The chair starts moving smoothly and quietly.



STEP 5

The chair automatically stops when it reaches the end of its journey at the top or bottom of the stairs.



STEP 6

Turn the seat very simply and dismount safely.



→ CONVENIENT

The unique controls are light and easy to operate.



→ SAFE

The safety system will stop the lift automatically if it encounters any obstacles.



→ CONVENIENT

The seat has a patented slide function to be as compact as possible.



THE CLEVER SOLUTIONS OF THE TWO

The TWO offers a good range of clever solutions. The seat can be swivelled either manually or automatically (depending on model) at the top of the stairs, allowing a safe dismount.

The TWO comes with an assisted manual folding footrest as standard but a motorised version is available for added convenience.



The TWO rail can be fully adapted to the interior of your home. Is there little space in front of your stairs? In that case you may opt for a sliding or foldable rail extension. Of course, your TWO can also be fitted with parking bends at the top and bottom of your stairs.




FOLDING EXTENSION
→ Folded in



→ Folded out



SLIDING EXTENSION
→ Slid in



→ Slid out

OTOLIFT **TWO**

Choose your upholstery
and paint colour



Upholstery



Grey
Fabric



Cream
Fabric



Red
Fabric



Bordeaux red
Fabric



Brown
Fabric



Violet
Fabric




Marine blue
Fabric



Pine green
Fabric



Variant beige
Fabric



Variant brown
Fabric




Hazelnut
Synthetic Leather



Reddish Brown
Synthetic Leather



Blue
Synthetic Leather



Black
Synthetic Leather



Taupe
Leather




Chestnut Brown
Leather


Paintwork



RAL 9001 Creamy white



RAL 7035 Light grey



RAL 8025 Pale brown





OTOLIFT TWO

Specials



CHILDREN

When it comes to stairlifts for children one issue is paramount: their safety. At Otolift you can rest assured that the stairlifts for children comply with the strictest requirements. There is an ideal Otolift for every age group, all based on one of our models. The stairlifts can be provided with separate controls for the parents, special seatbelts, electrically height adjustable footrests, and of course special seat elements.



INCREASED LIFTING CAPACITY - 165 KG

You may have gained weight as a result of using medicines over the past few years. Or perhaps as a result of reduced physical activity. Whatever the reason, you may need a stairlift built with more weight capacity.

Otolift will have the solution: we have developed a model suitable for up to 165 kg. This stairlift offers the same advantages as the standard version of the TWO, apart from its ability to carry 165 kg. What makes this lift so unique is the fact that the most important dimensions are the same. It is only the seat that is wider.

Every Otolift is build specifically to the users requirements. A number of alternative options are available to help ensure the lift completely meets your needs.

Whether it is for your child or for yourself, the Otolift Specials provide an excellent and above all safe solution!



WEATHER RESISTANT

Would you also like to experience the convenience of an Otolift outdoors? Once again no problem. There is an Otolift for straight or curved stairs. In the case of curved stairs only the rail is weather resistant; you need a dry spot for the chair. Thanks to the compact design these outdoor lifts cause no inconvenience to your neighbours. The level of comfort, safety and stability is just as high as that of Otolift stairlifts for indoor use.

Depicted is a Parallel stairlift, the Otolift straight unit.



STAND & PERCH CHAIR

Do you have knee problems which prevent you from fully bending your legs? Otolift offers a clever solution: a lift that can be used partially standing up so you hardly have to bend your knees at all. You 'lean' as it were on your legs, while sitting on a raised seat. This lift is also ideal for tall persons who have a narrow staircase in their home.

OTOLIFT TWO

Installation



← Quick and clean installation, in half a day

↓ Fixed to the staircase




About Otolift

Otolift is a Dutch family company, established in 1891. Since 1968, Otolift has been 100% producer of top quality stairlifts, and takes a lot of pride in being the benchmark in the industry.



→ 100% Produced by Otolift



→ 100% Designed by Otolift

→ Since 1891

Your Otolift approved dealer:

